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Bibliografia

1. Ilie N. (2018). Maturation of restorative glass ionomers with simplified application procedure. *Journal of dentistry*, 79, 46–52. <https://doi.org/10.1016/j.jdent.2018.09.008>
2. Ngo, H. C., Mount, G., Mc Intyre, J., Tuisuva, J., & Von Doussa, R. J. (2006). Chemical exchange between glass-ionomer restorations and residual carious dentine in permanent molars: an in vivo study. *Journal of dentistry*, 34(8), 608–613. <https://doi.org/10.1016/j.jdent.2005.12.012>
3. Ana Flávia Bissoto Calvo, Fabiana Bucholdz Teixeira Alves, Tathiane Larissa Lenzi, Tamara Kerber Tedesco, Alessandra Reis, Alessandro Dourado Loguercio, Daniela Prócida Raggio. Glass ionomer cements bond stability in caries-affected primary dentine. *International Journal of Adhesion & Adhesives* 48 (2014) 183–187
4. Francisconi, L. F., Scaffa, P. M., de Barros, V. R., Coutinho, M., & Francisconi, P. A. (2009). Glass ionomer cements and their role in the restoration of non-cariou cervical lesions. *Journal of applied oral science : revista FOB*, 17(5), 364–369. <https://doi.org/10.1590/s1678-775720090005000003>
5. Kim, Y. G., & Hirano, S. (1999). Setting shrinkage and hygroscopic expansion of resin-modified glass-ionomer in experimental cylindrical cavities. *Dental materials journal*, 18(1), 63–75. <https://doi.org/10.4012/dmj.18.63>
6. Cheetham, J. J., Palamara, J. E., Tyas, M. J., & Burrow, M. F. (2014). A comparison of resin-modified glass-ionomer and resin composite polymerisation shrinkage stress in a wet environment. *Journal of the mechanical behavior of biomedical materials*, 29, 33–41. <https://doi.org/10.1016/j.jmbbm.2013.07.003>
7. Naoum, S. J., Mutzelburg, P. R., Shumack, T. G., Thode, D., Martin, F. E., & Ellakwa, A. E. (2015). Reducing composite restoration polymerization shrinkage stress through resin modified glass-ionomer based adhesives. *Australian dental journal*, 60(4), 490–496. <https://doi.org/10.1111/adj.12265>
8. Perera, D., Yu, S., Zeng, H., Meyers, I. A., & Walsh, L. J. (2020). Acid Resistance of Glass Ionomer Cement Restorative Materials. *Bioengineering (Basel, Switzerland)*, 7(4), 150. <https://doi.org/10.3390/bioengineering7040150>
9. Tiwari, S., Kenchappa, M., Bhayya, D., Gupta, S., Saxena, S., Satyarth, S., Singh, A., & Gupta, M. (2016). Antibacterial Activity and Fluoride Release of Glass-Ionomer Cement, Compomer and Zirconia Reinforced Glass-Ionomer Cement. *Journal of clinical and diagnostic research : JCDR*, 10(4), ZC90–ZC93. <https://doi.org/10.7860/JCDR/2016/16282.7676>
10. Hafshejani, T. M., Zamanian, A., Venugopal, J. R., Rezvani, Z., Sefat, F., Saeb, M. R., Vahabi, H., Zarrintaj, P., & Mozafari, M. (2017). Antibacterial glass-ionomer cement restorative materials: A critical review on the current status of extended release formulations. *Journal of controlled release: official journal of the Controlled Release Society*, 262, 317–328. <https://doi.org/10.1016/j.jconrel.2017.07.041>
11. Cosgun, A., Bolgul, B., & Duran, N. (2019). In vitro investigation of antimicrobial effects, nanohardness, and cytotoxicity of different glass ionomer restorative materials in dentistry. *Nigerian journal of clinical practice*, 22(3), 422–431. https://doi.org/10.4103/njcp.njcp_429_18
12. Hoshika, S., De Munck, J., Sano, H., Sidhu, S. K., & Van Meerbeek, B. (2015). Effect of Conditioning and Aging on the Bond Strength and Interfacial Morphology of Glass-ionomer Cement Bonded to Dentine. *The journal of adhesive dentistry*, 17(2), 141–146. <https://doi.org/10.3290/j.jad.a33994>
13. Hoshika, S., Ting, S., Ahmed, Z., Chen, F., Toida, Y., Sakaguchi, N., Van Meerbeek, B., Sano, H., & Sidhu, S. K. (2021). Effect of conditioning and 1 year aging on the bond strength and interfacial morphology of glass-ionomer cement bonded to dentine. *Dental materials : official publication of the Academy of Dental Materials*, 37(1), 106–112. <https://doi.org/10.1016/j.dental.2020.10.016>
14. Hatirli, H., Yasa, B., & Çelik, E. U. (2021). Clinical performance of high-viscosity glass ionomer and resin composite on minimally invasive occlusal restorations performed without rubber-dam isolation: a two-year randomised split-mouth study. *Clinical oral investigations*, 25(9), 5493–5503. <https://doi.org/10.1007/s00784-021-03857-0>
15. Jhamak Nourmohammadi, Reza Salarian, Mehran Solati-Hashjin, Fatollah Moztarzadeh. Dissolution behavior and fluoride release from new glass composition used in glass ionomer cements. *Ceramics International* 33 (2007) 557–561
16. Schlafer, S., Bornmann, T., Paris, S., & Göstemeyer, G. (2021). The impact of glass ionomer cement and composite resin on microscale pH in cariogenic biofilms and demineralization of dental tissues. *Dental materials : official publication of the Academy of Dental Materials*, 37(10), 1576–1583. <https://doi.org/10.1016/j.dental.2021.08.007>
17. Featherstone, J. D., Glena, R., Shariati, M., & Shields, C. P. (1990). Dependence of in vitro demineralization of apatite and remineralization of dental enamel on fluoride concentration. *Journal of dental research*, 69 Spec No, 620–636. <https://doi.org/10.1177/00220345900690S121>
18. Gurgan, S., Kutuk, Z. B., Yalcin Cakir, F., & Ergin, E. (2020). A randomized controlled 10 years follow up of a glass ionomer restorative material in class I and class II cavities. *Journal of dentistry*, 94, 103175. <https://doi.org/10.1016/j.jdent.2019.07.013>
19. Miletić, I., Baraba, A., Basso, M., Pulcini, M. G., Marković, D., Perić, T., Ozkaya, C. A., & Turkun, L. S. (2020). Clinical Performance of a Glass-Hybrid System Compared with a Resin Composite in the Posterior Region: Results of a 2-year Multicenter Study. *The journal of adhesive dentistry*, 22(3), 235–247. <https://doi.org/10.3290/j.jad.a44547>